

The PPE on ProLens Frame Collection

Along with wearing a face mask, sealed safety eyewear can shield your patients from germs, airborne droplets and other harmful elements. This is the PPE on ProLens collection – sealed safety eyewear that's just \$39.95 when ordered on ProLens. Order your complimentary demo kit today.

Full Seal Frame Options



6001 + 6001C Armourx

Size: 62-16

Colors: Black, Bronze, Crystal, Grey
Must order 6001 frame and 6001C full seal cushion.



6002 + 6002C Armourx

Size: 53-18

Colors: Crystal, Grey
Must order 6002 frame and 6002C full seal cushion.



6005 Armourx

Size: 62-16

Colors: Black



6006 Armourx

Size: 60-17

Colors: Black, Blue, Grey



6008FS Armourx

Size: 59-19

Colors: Black



6009 Armourx

Size: 60-16

Colors: Black, Crystal, Grey



Baxter Rx Bolle Safety

Size: 63-18

Colors: Blue/Grey



Twister Bolle Safety

Size: 52-26

Colors: Black



WF590 Artcraft

Size: 57-19

Colors: Black/Blue, Violet/Grey, Black/Grey, Black Yellow, Navy/Orange

Full Seal Goggle Options



6007 Armourx

Size: 64-17

Colors: Black, Grey



6010 Armourx

Size: 64-17

Colors: Black, Brown



OG800 Onguard Safety

Size: 49-17

Colors: Black

To order your sample kit, contact your Walman Optical account manager!



WALMAN
OPTICAL

877.863.2759 | www.walmanoptical.com

Introducing: Fogless

Fogless is an anti-fog coating for any patient bothered by their glasses fogging up when wearing a face mask.

Facts about Fogless:

- A permanent, anti-fog coating applied as a dip hardcoat.
- No additional activators are needed! Fogless simply absorbs moisture so once the coating “fills up” it will tear the moisture out.
- Not compatible with non-glare.
- Can be applied to any lens material.
- Processed in-house to provide a profitable product for you with the quickest turntime possible.
- Included on ProLens at no additional charge and can be ordered with any ProLens job.
-

Maximize your Profitability: Use ProLens to Create Your Own PPE Package

Take ProLens one step further and use it to create a PPE package! Simply use ProLens to bundle a sealed safety frame with any lens and Fogless coating for the ultimate protective eyewear for your patients.



Order one of the sealed safety eyewear styles from the PPE on ProLens collection (just +\$39.95 on ProLens!)



Select any premium lens available on ProLens.



Choose Fogless anti-fog coating at no additional charge.



You'll have happy, healthy and protected patients.

**Contact your Walman Optical account manager for more
PPE on ProLens details.**